

Negatoscope Mammographe SEN'X Mammography viewer SEN'X

Modèle Ella **SEN'X SX8**

Caractéristiques techniques :

- nombre d'allumage : 1,
- surface utile : 73 x 43 cm,
- dimensions ext : 103 x 72 x 13 cm,
- poids 33 kg,
- nombre de tubes : 12 x 18 W
- **nombre de clichés :**

8 en 18 x 24 cm,

4 horizontales en 24 x 30 cm et 2 en 36 x 43 cm,

- luminance > 7000,
- écart d'homogénéité à 10 cm : 5%.

Characteristics

- Wall fixing: recessed holes for hooking at the back of the device.
- Light source: straight fluorescent tubes, "daylight" type, colour temperature > 6500 degrees Kelvin.
- Diffuser screen: acrylic, semi-translucent diffuser, thickness 3 mm.
- Lighting intensity: can be regulated from 1500 to 7000 cd / sq. m depending on model. Extremely uniform (standard DIN 6856).
- Film holder: treated steel rollers, non-wearing and non-loosening.
- Adjustable nylon wire to secure small formats.
- Central film grip (only for SEN'X 12)
- Lighting control: by electronic dimmer.
- Power supply cable: standard, 2 poles + earth.
- Mains electricity supply: 210/240V, 50/60Hz

Use:

Viewing area: 1 panel (73 x 43cm)

External dimension: 103 x 72 x 13cm

