

Cassette Capteur Plan Radiologie Numérique DR cassette Flat Panel Detector

FUJIFILM FDR D-EVO WIRELESS WiFi - G43i - 43*43cm - 17*17" - DR-ID602SE - S/N 37100313 - Année 2013 - Year 2013

+ Image Processing unit Console d'acquisition *FUJI Console Advance APL Software V7.1* - DR-ID300CL - PC HP Compaq Pro 6300 + DR-ID600MP / DR-ID600PU Power Supply Unit + Barco Touchscreen MDRC-1119TS
+ Interface Cable
+ D-Link Wireless access point AP - DAP 2553

Type	Cassette size detector with ISS (Irradiation Side Sampling system)
Scintillator	GOS (Gadolinium oxysulfide)
Detector external size	460 x 460 x 15 mm (Approx.) [18" x 18" x 0.6"]
Weight	4.0 kg [9lbs.] (including battery)
Pixel pitch	0.15 mm
Pixels	2816 x 2816 pixels
Wireless standard	IEEE 802.11n, 5.2GHz
Image preview	Approx. 1 sec
Cycle time	Approx. 8sec (wired mode) / Approx. 11 sec (wireless mode)
Battery recharging time	Approx. 3 hours
Battery performance	Standby: Approx. 3h 30min Number of exposures*: Approx. 500 exposures (@ 12sec cycle) <small>*When it's connected to the X-ray equipment directly</small>

TESTS VIDEO :

<https://youtu.be/LiA7JN7TZyg>



Ikonex Medical